

# Оборудование для культур клеток

## Технические характеристики

### По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231  
Ангарск (3955)60-70-56  
Архангельск (8182)63-90-72  
Астрахань (8512)99-46-04  
Барнаул (3852)73-04-60  
Белгород (4722)40-23-64  
Благовещенск (4162)22-76-07  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Владикавказ (8672)28-90-48  
Владимир (4922)49-43-18  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89  
Иваново (4932)77-34-06  
Ижевск (3412)26-03-58  
Иркутск (395)279-98-46  
Казань (843)206-01-48  
Россия +7(495)268-04-70

Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Коломна (4966)23-41-49  
Кострома (4942)77-07-48  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Курган (3522)50-90-47  
Липецк (4742)52-20-81  
Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41  
Нижний Новгород (831)429-08-12  
Новокузнецк (3843)20-46-81  
Ноябрьск (3496)41-32-12  
Новосибирск (383)227-86-73  
Киргизия +996(312)-96-26-47

Омск (3812)21-46-40  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Петрозаводск (8142)55-98-37  
Псков (8112)59-10-37  
Пермь (342)205-81-47  
Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Саранск (8342)22-96-24  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78  
Севастополь (8692)22-31-93  
Симферополь (3652)67-13-56  
Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Сургут (3462)77-98-35  
Казахстан +7(7172)727-132

Сыктывкар (8212)25-95-17  
Тамбов (4752)50-40-97  
Тверь (4822)63-31-35  
Тольятти (8482)63-91-07  
Томск (3822)98-41-53  
Тула (4872)33-79-87  
Тюмень (3452)66-21-18  
Ульяновск (8422)24-23-59  
Улан-Удэ (3012)59-97-51  
Уфа (347)229-48-12  
Хабаровск (4212)92-98-04  
Чебоксары (8352)28-53-07  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Чита (3022)38-34-83  
Якутск (4112)23-90-97  
Ярославль (4852)69-52-93

# Corning® Cap Inserts for Cryogenic Vials



- Cap inserts provide color coding for easy sample identification.
- Inserts are packaged in resealable bags.
- Nonsterile
- Cap inserts fit all Corning® cryogenic vials.

Color

|  |          |
|--|----------|
|  | Assorted |
|  | Blue     |
|  | Green    |
|  | Red      |
|  | White    |
|  | Yellow   |

Products

|                |                                                         |
|----------------|---------------------------------------------------------|
| Product Number | 2015                                                    |
| Product Name   | Corning® White Polypropylene Cryogenic Vial Cap Inserts |
| Qty./Pk        | 50 / Pk                                                 |
| Qty./Cs        | 500 / Cs                                                |

|                |                                                            |
|----------------|------------------------------------------------------------|
| Product Number | 2016                                                       |
| Product Name   | Corning® Blue Polypropylene Cryogenic Vial Cap Inserts     |
| Qty./Pk        | 50 / Pk                                                    |
| Qty./Cs        | 500 / Cs                                                   |
| Product Number | 2017                                                       |
| Product Name   | Corning® Red Polypropylene Cryogenic Vial Cap Inserts      |
| Qty./Pk        | 50 / Pk                                                    |
| Qty./Cs        | 500 / Cs                                                   |
| Product Number | 2018                                                       |
| Product Name   | Corning® Green Polypropylene Cryogenic Vial Cap Inserts    |
| Qty./Pk        | 50 / Pk                                                    |
| Qty./Cs        | 500 / Cs                                                   |
| Product Number | 2019                                                       |
| Product Name   | Corning® Yellow Polypropylene Cryogenic Vial Cap Inserts   |
| Qty./Pk        | 50 / Pk                                                    |
| Qty./Cs        | 500 / Cs                                                   |
| Product Number | 430499                                                     |
| Product Name   | Corning® Assorted Polypropylene Cryogenic Vial Cap Inserts |
| Qty./Pk        | 50 / Pk                                                    |
| Qty./Cs        | 500 / Cs                                                   |

# Corning® External Thread Cryogenic Vials



- Manufactured from polypropylene to withstand temperatures down to -196°C
- Black graduations with large marking spot
- Vials have a silicone washer for a secure seal.
- 2 mL self-standing vials are available in five cap colors
- Vials may be color coded with inserts.
- Self-standing vials have a special base design allowing them to be locked into cryogenic rack and tray (Cat. No. 431131) for single-handed manipulation.
- Free foam rack with each case
- RNase-/DNase-free
- Sterile
- Nonpyrogenic

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquefied gas. Always use appropriate safety equipment when removing vials from cryogenic storage.

## Vial Volume

|  |        |
|--|--------|
|  | 1.2 mL |
|  | 2 mL   |
|  | 4 mL   |
|  | 5 mL   |

## Bottom Style

|  |         |
|--|---------|
|  | Conical |
|  | Round   |

## Self-standing

|  |     |
|--|-----|
|  | No  |
|  | Yes |


## Products



|                       |                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------|
| <b>Product Number</b> | 430658                                                                                            |
| <b>Product Name</b>   | Corning® 1.2 mL External Threaded Polypropylene Cryogenic Vial, Self-Standing with Conical Bottom |
| <b>Qty./Pk</b>        | 50 / Pk                                                                                           |
| <b>Qty./Cs</b>        | 500 / Cs                                                                                          |
| <b>Product Number</b> | 430659                                                                                            |
| <b>Product Name</b>   | Corning® 2 mL External Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom     |
| <b>Qty./Pk</b>        | 50 / Pk                                                                                           |
| <b>Qty./Cs</b>        | 500 / Cs                                                                                          |
| <b>Product Number</b> | 430661                                                                                            |
| <b>Product Name</b>   | Corning® 2 mL External Threaded Polypropylene Cryogenic Vial with Round Bottom                    |
| <b>Qty./Pk</b>        | 50 / Pk                                                                                           |
| <b>Qty./Cs</b>        | 500 / Cs                                                                                          |
| <b>Product Number</b> | 430662                                                                                            |
| <b>Product Name</b>   | Corning® 4 mL External Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom     |
| <b>Qty./Pk</b>        | 50 / Pk                                                                                           |
| <b>Qty./Cs</b>        | 500 / Cs                                                                                          |
| <b>Product Number</b> | 430663                                                                                            |
| <b>Product Name</b>   | Corning® 5 mL External Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom     |
| <b>Qty./Pk</b>        | 50 / Pk                                                                                           |
| <b>Qty./Cs</b>        | 500 / Cs                                                                                          |

|                |                                                                                                                                 |
|----------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Number | 430674                                                                                                                          |
| Product Name   | Corning® 5 mL External Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom, Cap Unassembled, Unprinted, Bulk |
| Qty./Pk        | 125 / Pk                                                                                                                        |
| Qty./Cs        | 500 / Cs                                                                                                                        |

# Corning® Cryogenic Vial Racks and Storage Boxes



- Reusable racks are designed for use with most cryogenic vials
- Cat. No. 431131 has a locking feature for use with all Corning self-standing vials
- Cat. No. 431121 accepts internally threaded cryogenic vials only.

## Cryogenic Storage Format

|                                  |
|----------------------------------|
| 1-2 mL Internally-threaded Vials |
| 1-2 mL Vials                     |
| 4-5 mL Vials                     |

## Products

|                |                                                                            |
|----------------|----------------------------------------------------------------------------|
| Product Number | 431119                                                                     |
| Product Name   | Corning® Polycarbonate 1 - 2 mL Cryogenic Vial Storage Box, Holds 81 Vials |
| Qty./Pk        | 5 / Pk                                                                     |
| Qty./Cs        | 10 / Cs                                                                    |
| Product Number | 431120                                                                     |
| Product Name   | Corning® Polycarbonate 4 - 5 mL Cryogenic Vial Storage Box, Holds 81 Vials |
| Qty./Pk        | 5 / Pk                                                                     |
| Qty./Cs        | 10 / Cs                                                                    |

|                |                                                                             |
|----------------|-----------------------------------------------------------------------------|
| Product Number | 431121                                                                      |
| Product Name   | Corning® Polycarbonate 1 - 2 mL Cryogenic Vial Storage Box, Holds 100 Vials |
| Qty./Pk        | 5 / Pk                                                                      |
| Qty./Cs        | 10 / Cs                                                                     |

# Corning® Polypropylene Cryogenic Vial Rack, Holds 50 Vials

Product Number 431131



- Reusable polypropylene orange rack holds 50 cryogenic vials
- Has locking feature for use with all Corning® self-standing vials; allows one handed opening

## Details

|                |          |
|----------------|----------|
| Product Number | 431131   |
| Qty./Pk        | 2 / Pk   |
| Qty./Cs        | 2 / Cs   |
| Brand          | Corning® |

# Corning® Bar Coded External and Internal Threaded Polypropylene Cryogenic Vials



Threading

|          |
|----------|
| External |
| Internal |

Barcode

|           |
|-----------|
| 1D and 2D |
| 1D        |

Products

|                |                                                                                       |
|----------------|---------------------------------------------------------------------------------------|
| Product Number | 8670                                                                                  |
| Product Name   | Corning® 2 mL Internal Threaded Polypropylene Cryogenic Vial with 1D and 2D Bar codes |
| Qty./Pk        | 50 / Pk                                                                               |
| Qty./Cs        | 500 / Cs                                                                              |

|                       |                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------|
| <b>Product Number</b> | 8671                                                                                  |
| <b>Product Name</b>   | Corning® 2 mL External Threaded Polypropylene Cryogenic Vial with 1D and 2D Bar codes |
| <b>Qty./Pk</b>        | 50 / Pk                                                                               |
| <b>Qty./Cs</b>        | 500 / Cs                                                                              |
| <b>Product Number</b> | 8672                                                                                  |
| <b>Product Name</b>   | Corning® 2 mL Internal Threaded Polypropylene Cryogenic Vial with 1D Barcode          |
| <b>Qty./Pk</b>        | 50 / Pk                                                                               |
| <b>Qty./Cs</b>        | 500 / Cs                                                                              |
| <b>Product Number</b> | 8676                                                                                  |
| <b>Product Name</b>   | Corning® 2 mL External Threaded Polypropylene Cryogenic Vial with 1D Barcode          |
| <b>Qty./Pk</b>        | 50 / Pk                                                                               |
| <b>Qty./Cs</b>        | 500 / Cs                                                                              |

# Corning® Internal Thread Cryogenic Vials



- Manufactured from polypropylene to withstand temperatures down to -196°C
- Black graduations with large marking spot
- Vials have a silicone washer for a secure seal.
- 2 mL self-standing vials are available in five cap colors.
- Vials may be color coded with inserts.
- Self-standing vials have a special base design allowing them to be locked into cryogenic rack and tray (Cat. No. 431131) for single-handed manipulation.
- Free foam rack with each case
- RNase-/DNase-free
- Sterile
- Nonpyrogenic

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquefied gas. Always use appropriate safety equipment when removing vials from cryogenic storage.

## Color

|  |          |
|--|----------|
|  | Assorted |
|  | Blue     |
|  | Green    |
|  | Orange   |
|  | Red      |
|  | White    |
|  | Yellow   |

## Capacity

|  |        |
|--|--------|
|  | 1.2 mL |
|  | 2 mL   |
|  | 4 mL   |
|  | 5 mL   |

## Self-standing



|     |
|-----|
| No  |
| Yes |

Qty/Cs

|        |
|--------|
| 250/Cs |
| 500/Cs |

Products

|                |                                                                                                   |
|----------------|---------------------------------------------------------------------------------------------------|
| Product Number | 430487                                                                                            |
| Product Name   | Corning® 1.2 mL Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Conical Bottom |
| Qty./Pk        | 50 / Pk                                                                                           |
| Qty./Cs        | 500 / Cs                                                                                          |

|                |                                                                                               |
|----------------|-----------------------------------------------------------------------------------------------|
| Product Number | 430488                                                                                        |
| Product Name   | Corning® 2 mL Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom |
| Qty./Pk        | 50 / Pk                                                                                       |
| Qty./Cs        | 500 / Cs                                                                                      |

|                |                                                                                |
|----------------|--------------------------------------------------------------------------------|
| Product Number | 430489                                                                         |
| Product Name   | Corning® 2 mL Internal Threaded Polypropylene Cryogenic Vial with Round Bottom |
| Qty./Pk        | 50 / Pk                                                                        |
| Qty./Cs        | 500 / Cs                                                                       |

|                |                                                                                |
|----------------|--------------------------------------------------------------------------------|
| Product Number | 430490                                                                         |
| Product Name   | Corning® 4 mL Internal Threaded Polypropylene Cryogenic Vial with Round Bottom |
| Qty./Pk        | 50 / Pk                                                                        |
| Qty./Cs        | 500 / Cs                                                                       |

|                |                                                                                                               |
|----------------|---------------------------------------------------------------------------------------------------------------|
| Product Number | 430491                                                                                                        |
| Product Name   | Corning® 4 mL Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom                 |
| Qty./Pk        | 50 / Pk                                                                                                       |
| Qty./Cs        | 500 / Cs                                                                                                      |
| Product Number | 430492                                                                                                        |
| Product Name   | Corning® 5 mL Internal Threaded Polypropylene Cryogenic Vial with Round Bottom                                |
| Qty./Pk        | 50 / Pk                                                                                                       |
| Qty./Cs        | 500 / Cs                                                                                                      |
| Product Number | 430656                                                                                                        |
| Product Name   | Corning® 5 mL Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom                 |
| Qty./Pk        | 50 / Pk                                                                                                       |
| Qty./Cs        | 500 / Cs                                                                                                      |
| Product Number | 431386                                                                                                        |
| Product Name   | Corning® 2 mL Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom                 |
| Qty./Pk        | 50 / Pk                                                                                                       |
| Qty./Cs        | 250 / Cs                                                                                                      |
| Product Number | 431416                                                                                                        |
| Product Name   | Corning® 2 mL Assorted Colors Internal Threaded Polypropylene Cryogenic Vial, Self-Standing with Round Bottom |
| Qty./Pk        | 50 / Pk                                                                                                       |
| Qty./Cs        | 500 / Cs                                                                                                      |

|                |                                                                                                              |
|----------------|--------------------------------------------------------------------------------------------------------------|
| Product Number | 431417                                                                                                       |
| Product Name   | Corning® 2 mL Yellow Cap Internal Threaded Polypropylene Cryogenic Vial, Self-<br>Standing with Round Bottom |
| Qty./Pk        | 50 / Pk                                                                                                      |
| Qty./Cs        | 500 / Cs                                                                                                     |

# Corning® CoolCell® Containers



Corning CoolCell, in combination with a -80°C freezer, will provide alcohol-free freezing at the rate of -1°C/minute that is ideal for cryopreservation of most cells and cell lines. Using a combination of uniform-density cross-linked polyethylene foam, a solid state core, and radial vial symmetry, freezing profiles are consistent and reproducible. The foam is non-absorbent and will impose negligible change in the freezer environment; thereby protecting nearby frozen samples. The removal of vials is accomplished by tipping the device upside down so the vials slide out.

Features: Unique features of controlled-rate freezing with Corning® CoolCell® include:

- Ease of use
- Alcohol and fluid-free freezing
- Lower cost of use than alcohol-based devices
- High cell recovery and cell viability
- Reproducibility
- Simple, consistent way to standardize controlled-rate freezing

## Cryogenic Storage Format

|                                 |
|---------------------------------|
| 12 x 1 or 2 mL Cryogenic Vials  |
| 12 x 3.5 - 5 mL Cryogenic Vials |
| 30 x 1 or 2 mL Cryogenic Vials  |
| 12 x 2 mL Serum Vials           |
| 6 x 10 mL Serum Vials           |

## Color

|                             |
|-----------------------------|
| Green                       |
| Orange                      |
| Pink                        |
| Purple                      |
| purple, orange, green, pink |

## Products

|                |                                                                                               |
|----------------|-----------------------------------------------------------------------------------------------|
| Product Number | 432000                                                                                        |
| Product Name   | Corning® CoolCell® Cell Freezing Container, for 12 x 1 mL or 2 mL Cryogenic Vials, Purple     |
| Qty./Pk        | 1 / Pk                                                                                        |
| Qty./Cs        | 1 / Cs                                                                                        |
| Product Number | 432001                                                                                        |
| Product Name   | Corning® CoolCell® LX, Cell Freezing Container, for 12 x 1 mL or 2 mL Cryogenic Vials, Purple |
| Qty./Pk        | 1 / Pk                                                                                        |
| Qty./Cs        | 1 / Cs                                                                                        |
| Product Number | 432002                                                                                        |
| Product Name   | Corning® CoolCell® LX, Cell Freezing Container, for 12 x 1 mL or 2 mL Cryogenic Vials, Green  |
| Qty./Pk        | 1 / Pk                                                                                        |
| Qty./Cs        | 1 / Cs                                                                                        |
| Product Number | 432003                                                                                        |
| Product Name   | Corning® CoolCell® LX, Cell Freezing Container, for 12 x 1 mL or 2 mL Cryogenic Vials, Orange |
| Qty./Pk        | 1 / Pk                                                                                        |
| Qty./Cs        | 1 / Cs                                                                                        |
| Product Number | 432004                                                                                        |
| Product Name   | Corning® CoolCell® LX, Cell Freezing Container, for 12 x 1 mL or 2 mL Cryogenic Vials, Pink   |
| Qty./Pk        | 1 / Pk                                                                                        |
| Qty./Cs        | 1 / Cs                                                                                        |

|                       |                                                                                                      |
|-----------------------|------------------------------------------------------------------------------------------------------|
| <b>Product Number</b> | 432005                                                                                               |
| <b>Product Name</b>   | Corning® CoolCell® 5 mL LX, Cell Freezing Container, for 12 x 3.5 mL to 5 mL Cryogenic Vials, Purple |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                               |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                               |
| <b>Product Number</b> | 432006                                                                                               |
| <b>Product Name</b>   | Corning® CoolCell® FTS30, Freezing Container, for 30 x 1 mL or 2 mL Cryogenic Vials, Purple          |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                               |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                               |
| <b>Product Number</b> | 432007                                                                                               |
| <b>Product Name</b>   | Corning® CoolCell® FTS30 Freezing Container, for 30 x 1 mL or 2 mL Cryogenic Vials, Orange           |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                               |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                               |
| <b>Product Number</b> | 432008                                                                                               |
| <b>Product Name</b>   | Corning® CoolCell® FTS30, Freezing Container, for 30 x 1 mL or 2 mL Cryogenic Vials, Green           |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                               |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                               |
| <b>Product Number</b> | 432009                                                                                               |
| <b>Product Name</b>   | Corning® CoolCell® FTS30, Freezing Container, for 30 x 1 mL or 2 mL Cryogenic Vials, Pink            |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                               |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                               |

# Corning® Cryogenic Storage Boxes for Bar Coded Cryogenic Vials

Vial Number

|    |
|----|
| 81 |
|----|

|     |
|-----|
| 100 |
|-----|

## Products

|                |                                                                                         |
|----------------|-----------------------------------------------------------------------------------------|
| Product Number | 8673                                                                                    |
| Product Name   | Corning® Polycarbonate 1 - 2 mL 2D barcoded Cryogenic Vial Storage Box, Holds 81 Vials  |
| Qty./Pk        | 5 / Pk                                                                                  |
| Qty./Cs        | 10 / Cs                                                                                 |
| Product Number | 8674                                                                                    |
| Product Name   | Corning® Polycarbonate 1 - 2 mL 2D barcoded Cryogenic Vial Storage Box, Holds 100 Vials |
| Qty./Pk        | 5 / Pk                                                                                  |
| Qty./Cs        | 10 / Cs                                                                                 |

# Corning® Cryogenic Vial Grippers, Case of 5, Assorted Colors

Product Number 432136



Cryogenic vial grippers feature a unique design to grasp internal or external-thread cryogenic vials. Grippers enable to easily sort or move vials while maintaining sterility and protecting fingers from frozen vials, dry ice and liquid nitrogen. 5 per case.

## Details

|                |             |
|----------------|-------------|
| Product Number | 432136      |
| Qty./Pk        | 1 / Pk      |
| Qty./Cs        | 5 / Cs      |
| Brand          | Corning®    |
| Application    | Diagnostics |



# Corning® CoolCell® Filler Vials



Corning CoolCell, in combination with a -80°C freezer, will provide alcohol-free freezing at the rate of -1°C/minute that is ideal for cryopreservation of most cells and cell lines. Using a combination of uniform-density cross-linked polyethylene foam, a solid state core, and radial vial symmetry, freezing profiles are consistent and reproducible. The foam is non-absorbent and will impose negligible change in the freezer environment; thereby protecting nearby frozen samples. The removal of vials is accomplished by tipping the device upside down so the vials slide out.

Features: Unique features of controlled-rate freezing with Corning® CoolCell® include:

- Ease of use
- Alcohol and fluid-free freezing
- Lower cost of use than alcohol-based devices
- High cell recovery and cell viability
- Reproducibility
- Simple, consistent way to standardize controlled-rate freezing

**Vial Volume**

|  |      |
|--|------|
|  | 2 mL |
|  | 5 mL |


## Products

|                |                                                                                             |
|----------------|---------------------------------------------------------------------------------------------|
| Product Number | 432076                                                                                      |
| Product Name   | Corning® CoolCell® 2 mL Filler Vials for Use with CoolCell LX and CoolCell FTS30 Containers |
| Qty./Pk        | 1 / Pk                                                                                      |
| Qty./Cs        | 6 / Cs                                                                                      |

|                |                                                                              |
|----------------|------------------------------------------------------------------------------|
| Product Number | 432077                                                                       |
| Product Name   | Corning® CoolCell® 5 mL Filler Vials for Use with CoolCell 5 mL LX Container |
| Qty./Pk        | 1 / Pk                                                                       |
| Qty./Cs        | 6 / Cs                                                                       |

# Removable Cryogenic vial Module for Corning® CoolCell® FTS30

Product Number 432078



Corning CoolCell FTS30 Vial Module is a holder for 30 1.0 mL or 2.0 mL cryogenic vials that allows one-step insertion and removal of all 30 vials at once. It fits into a standard 5.0 x 5.0 x 2.0 inch cryostorage box and is compatible with dry ice and liquid nitrogen.

## Details

|                |          |
|----------------|----------|
| Product Number | 432078   |
| Qty./Pk        | 5 / Pk   |
| Qty./Cs        | 10 / Cs  |
| Brand          | Corning® |

# Corning® CoolRack® M Microcentrifuge Tube Modules



Color

|  |        |
|--|--------|
|  | Green  |
|  | Grey   |
|  | Orange |

Cryogenic Storage Format

|  |                                                                                        |
|--|----------------------------------------------------------------------------------------|
|  | 6 x 1.5 or 2 mL Microcentrifuge Tubes                                                  |
|  | 15 x 1.5 or 2 mL Microcentrifuge Tubes                                                 |
|  | 30 x 1.5 or 2 mL Microcentrifuge Tubes                                                 |
|  | 90 x 1.5 or 2 mL Microcentrifuge Tubes                                                 |
|  | 96 x 1.5 or 2 mL Microcentrifuge Tubes with A to H and 1 to 12 Row and Column Indexing |

Product Type

|  |                 |
|--|-----------------|
|  | CoolRack M15    |
|  | CoolRack M30    |
|  | CoolRack M6     |
|  | CoolRack M90    |
|  | CoolRack M96 ID |

Products

|                |                                                                            |
|----------------|----------------------------------------------------------------------------|
| Product Number | 432034                                                                     |
| Product Name   | Corning® CoolRack M6, Holds 6 x 1.5 or 2 mL Microcentrifuge Tubes, Gray    |
| Qty./Pk        | 1 / Pk                                                                     |
| Qty./Cs        | 1 / Cs                                                                     |
| Product Number | 432035                                                                     |
| Product Name   | Corning® CoolRack M6, Holds 6 x 1.5 or 2 mL Microcentrifuge Tubes, Green   |
| Qty./Pk        | 1 / Pk                                                                     |
| Qty./Cs        | 1 / Cs                                                                     |
| Product Number | 432036                                                                     |
| Product Name   | Corning® CoolRack M6, Holds 6 x 1.5 or 2 mL Microcentrifuge Tubes, Orange  |
| Qty./Pk        | 1 / Pk                                                                     |
| Qty./Cs        | 1 / Cs                                                                     |
| Product Number | 432037                                                                     |
| Product Name   | Corning® CoolRack M15, Holds 15 x 1.5 or 2 mL Microcentrifuge Tubes, Gray  |
| Qty./Pk        | 1 / Pk                                                                     |
| Qty./Cs        | 1 / Cs                                                                     |
| Product Number | 432038                                                                     |
| Product Name   | Corning® CoolRack M15, Holds 15 x 1.5 or 2 mL Microcentrifuge Tubes, Green |
| Qty./Pk        | 1 / Pk                                                                     |
| Qty./Cs        | 1 / Cs                                                                     |

|                       |                                                                             |
|-----------------------|-----------------------------------------------------------------------------|
| <b>Product Number</b> | 432039                                                                      |
| <b>Product Name</b>   | Corning® CoolRack M15, Holds 15 x 1.5 or 2 mL Microcentrifuge Tubes, Orange |
| <b>Qty./Pk</b>        | 1 / Pk                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                      |
| <b>Product Number</b> | 432041                                                                      |
| <b>Product Name</b>   | Corning® CoolRack M30, Holds 30 x 1.5 or 2 mL Microcentrifuge Tubes, Gray   |
| <b>Qty./Pk</b>        | 1 / Pk                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                      |
| <b>Product Number</b> | 432042                                                                      |
| <b>Product Name</b>   | Corning® CoolRack M30, Holds 30 x 1.5 or 2 mL Microcentrifuge Tubes, Green  |
| <b>Qty./Pk</b>        | 1 / Pk                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                      |
| <b>Product Number</b> | 432043                                                                      |
| <b>Product Name</b>   | Corning® CoolRack M30, Holds 30 x 1.5 or 2 mL Microcentrifuge Tubes, Orange |
| <b>Qty./Pk</b>        | 1 / Pk                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                      |
| <b>Product Number</b> | 432044                                                                      |
| <b>Product Name</b>   | Corning® CoolRack M90, Holds 90 x 1.5 or 2 mL Microcentrifuge Tubes, Gray   |
| <b>Qty./Pk</b>        | 1 / Pk                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                      |

# Corning® Ice Pans, Rectangular



Color

|  |        |
|--|--------|
|  | Blue   |
|  | Green  |
|  | Orange |
|  | Pink   |
|  | Purple |

Capacity

|  |     |
|--|-----|
|  | 1 L |
|  | 4 L |
|  | 9 L |

Product Type

|  |                        |
|--|------------------------|
|  | Ice pan, maxi          |
|  | Ice pan, midi          |
|  | Ice pan with lid, midi |
|  | Ice pan, mini          |

---

## Products

|                |                                                    |
|----------------|----------------------------------------------------|
| Product Number | 432093                                             |
| Product Name   | Corning® Ice Pan, Rectangular, Maxi 9L, Green      |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 1 / Cs                                             |
| Product Number | 432094                                             |
| Product Name   | Corning® Ice Pan, Rectangular, Maxi 9L, Blue       |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 1 / Cs                                             |
| Product Number | 432095                                             |
| Product Name   | Corning® Ice Pan, Rectangular, Maxi 9L, Red        |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 1 / Cs                                             |
| Product Number | 432096                                             |
| Product Name   | Corning® Ice Pan, Rectangular, Maxi 9L, Orange     |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 1 / Cs                                             |
| Product Number | 432097                                             |
| Product Name   | Corning® Ice Pan, Rectangular, Maxi 9L, Lime Green |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 1 / Cs                                             |
| Product Number | 432098                                             |
| Product Name   | Corning® Ice Pan, Rectangular, Maxi 9L, Pink       |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 1 / Cs                                             |



|                       |                                                             |
|-----------------------|-------------------------------------------------------------|
| <b>Product Number</b> | 432099                                                      |
| <b>Product Name</b>   | Corning® Ice Pan, Rectangular, Maxi 9L, Purple              |
| <b>Qty./Pk</b>        | 1 / Pk                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                      |
| <b>Product Number</b> | 432100                                                      |
| <b>Product Name</b>   | Corning® Ice Pan, Rectangular with Lid, Maxi 9L, Blue       |
| <b>Qty./Pk</b>        | 1 / Pk                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                      |
| <b>Product Number</b> | 432101                                                      |
| <b>Product Name</b>   | Corning® Ice Pan, Rectangular with Lid, Maxi 9L, Lime Green |
| <b>Qty./Pk</b>        | 1 / Pk                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                      |
| <b>Product Number</b> | 432102                                                      |
| <b>Product Name</b>   | Corning® Ice Pan, Rectangular with Lid, Maxi 9L, Purple     |
| <b>Qty./Pk</b>        | 1 / Pk                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                      |

# Corning® CoolBox™ Accessories



Color

|  |        |
|--|--------|
|  | Black  |
|  | Blue   |
|  | Green  |
|  | Grey   |
|  | Orange |

Product Type

|  |                                          |
|--|------------------------------------------|
|  | CoolRack Elastic Sleeves                 |
|  | CoolRack Temperature Strip, 1 °C to 8 °C |
|  | Extension collar for CoolBox 2XT         |
|  | Extension collar for CoolBox XT          |
|  | XT Cooling Core                          |

Products

|                |                                                                        |
|----------------|------------------------------------------------------------------------|
| Product Number | 432081                                                                 |
| Product Name   | XT Cooling Core for Corning® CoolBox™ XT Workstations, for 0.5° to 4°C |
| Qty./Pk        | 1 / Pk                                                                 |
| Qty./Cs        | 1 / Cs                                                                 |

|                |                                                                       |
|----------------|-----------------------------------------------------------------------|
| Product Number | 432082                                                                |
| Product Name   | XT Freezing Core for Corning® CoolBox™ XT Workstations, for Below 0°C |
| Qty./Pk        | 1 / Pk                                                                |
| Qty./Cs        | 1 / Cs                                                                |
| Product Number | 432083                                                                |
| Product Name   | Extension Collar for Corning® CoolBox™ XT, Purple                     |
| Qty./Pk        | 1 / Pk                                                                |
| Qty./Cs        | 1 / Cs                                                                |
| Product Number | 432084                                                                |
| Product Name   | Extension Collar for Corning® CoolBox™ XT, Green                      |
| Qty./Pk        | 1 / Pk                                                                |
| Qty./Cs        | 1 / Cs                                                                |
| Product Number | 432085                                                                |
| Product Name   | Extension Collar for Corning® CoolBox™ XT, Orange                     |
| Qty./Pk        | 1 / Pk                                                                |
| Qty./Cs        | 1 / Cs                                                                |
| Product Number | 432086                                                                |
| Product Name   | Extension Collar for Corning® CoolBox™ XT, Pink                       |
| Qty./Pk        | 1 / Pk                                                                |
| Qty./Cs        | 1 / Cs                                                                |
| Product Number | 432087                                                                |
| Product Name   | Extension Collar for Corning® CoolBox™ 2XT, Purple                    |
| Qty./Pk        | 1 / Pk                                                                |
| Qty./Cs        | 1 / Cs                                                                |

|                       |                                                    |
|-----------------------|----------------------------------------------------|
| <b>Product Number</b> | 432088                                             |
| <b>Product Name</b>   | Extension Collar for Corning® CoolBox™ 2XT, Green  |
| <b>Qty./Pk</b>        | 1 / Pk                                             |
| <b>Qty./Cs</b>        | 1 / Cs                                             |
| <b>Product Number</b> | 432089                                             |
| <b>Product Name</b>   | Extension Collar for Corning® CoolBox™ 2XT, Orange |
| <b>Qty./Pk</b>        | 1 / Pk                                             |
| <b>Qty./Cs</b>        | 1 / Cs                                             |
| <b>Product Number</b> | 432090                                             |
| <b>Product Name</b>   | Extension Collar for Corning® CoolBox™ 2XT, Pink   |
| <b>Qty./Pk</b>        | 1 / Pk                                             |
| <b>Qty./Cs</b>        | 1 / Cs                                             |

# Corning® CoolRack® M-PF Conical Tube Modules



**Product Type**

|                        |
|------------------------|
| CoolRack 500 µL M30-PF |
| CoolRack M15-PF        |
| CoolRack M30-PF        |

**Cryogenic Storage Format**

|                           |
|---------------------------|
| 15 x 1.5 mL conical tubes |
| 30 x 1.5 mL conical tubes |
| 30 x 500 µL conical tubes |

**Products**

|                |                                                                                                               |
|----------------|---------------------------------------------------------------------------------------------------------------|
| Product Number | 432046                                                                                                        |
| Product Name   | Corning® CoolRack 500ul M30-PF, Holds 30 x 500uL Microcentrifuge Tubes, Tapered Wells for Conical Tubes, Gray |
| Qty./Pk        | 1 / Pk                                                                                                        |
| Qty./Cs        | 1 / Cs                                                                                                        |

|                       |                                                                                                         |
|-----------------------|---------------------------------------------------------------------------------------------------------|
| <b>Product Number</b> | 432047                                                                                                  |
| <b>Product Name</b>   | Corning® CoolRack M15-PF, Holds 15 x 1.5mL Microcentrifuge Tubes, Tapered Wells for Conical Tubes, Gray |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                                  |
| <b>Product Number</b> | 432048                                                                                                  |
| <b>Product Name</b>   | Corning® CoolRack M30-PF, Holds 30 x 1.5mL Microcentrifuge Tubes, Tapered Wells for Conical Tubes, Gray |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                                  |

# Corning® CoolRack® 96 Storage Tube Modules



Product Type

|  |                      |
|--|----------------------|
|  | CoolRack 96 x 0.5 mL |
|  | CoolRack 96 x 1 mL   |

Cryogenic Storage Format

|  |                               |
|--|-------------------------------|
|  | 96 x 0.5 mL 2D bar code tubes |
|  | 96 x 1.4 mL 2D bar code tubes |

Products

|                |                                                       |
|----------------|-------------------------------------------------------|
| Product Number | 432056                                                |
| Product Name   | Corning® CoolRack 96x0.5mL, Holds 96 x 0.5mL 2D Tubes |
| Qty./Pk        | 1 / Pk                                                |
| Qty./Cs        | 1 / Cs                                                |
| Product Number | 432057                                                |
| Product Name   | Corning® CoolRack 96x1 mL, Holds 96 x 1.4 mL 2D Tubes |
| Qty./Pk        | 1 / Pk                                                |
| Qty./Cs        | 1 / Cs                                                |

# Corning® CoolBox™ XT Modules



Product Type

CoolBox 2XT

CoolBox XT

Cryogenic Storage Format

Holder and one Cooling Core

Holder and two Cooling Cores

Color

Green

Orange

Pink

Purple

## Products

|                |                                                                                        |
|----------------|----------------------------------------------------------------------------------------|
| Product Number | 432021                                                                                 |
| Product Name   | Corning® CoolBox™ XT All-day Cooling and Freezing Workstation, Single Capacity, Purple |
| Qty./Pk        | 1 / Pk                                                                                 |
| Qty./Cs        | 1 / Cs                                                                                 |



|                       |                                                                                         |
|-----------------------|-----------------------------------------------------------------------------------------|
| <b>Product Number</b> | 432022                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ XT All-day Cooling and Freezing Workstation, Single Capacity, Green   |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432023                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ XT All-day Cooling and Freezing Workstation, Single Capacity, Orange  |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432024                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ XT, Pink All-day Cooling and Freezing Workstation, Single Capacity    |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432025                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Purple |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432026                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Green  |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |

|                |                                                                                         |
|----------------|-----------------------------------------------------------------------------------------|
| Product Number | 432027                                                                                  |
| Product Name   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Orange |
| Qty./Pk        | 1 / Pk                                                                                  |
| Qty./Cs        | 1 / Cs                                                                                  |
| Product Number | 432028                                                                                  |
| Product Name   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Pink   |
| Qty./Pk        | 1 / Pk                                                                                  |
| Qty./Cs        | 1 / Cs                                                                                  |

# Corning® CoolBox™ XT Modules



Product Type

CoolBox 2XT

CoolBox XT

Cryogenic Storage Format

Holder and one Cooling Core

Holder and two Cooling Cores

Color

Green

Orange

Pink

Purple

## Products

|                |                                                                                        |
|----------------|----------------------------------------------------------------------------------------|
| Product Number | 432021                                                                                 |
| Product Name   | Corning® CoolBox™ XT All-day Cooling and Freezing Workstation, Single Capacity, Purple |
| Qty./Pk        | 1 / Pk                                                                                 |
| Qty./Cs        | 1 / Cs                                                                                 |

|                       |                                                                                         |
|-----------------------|-----------------------------------------------------------------------------------------|
| <b>Product Number</b> | 432022                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ XT All-day Cooling and Freezing Workstation, Single Capacity, Green   |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432023                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ XT All-day Cooling and Freezing Workstation, Single Capacity, Orange  |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432024                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ XT, Pink All-day Cooling and Freezing Workstation, Single Capacity    |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432025                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Purple |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |
| <b>Product Number</b> | 432026                                                                                  |
| <b>Product Name</b>   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Green  |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                  |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                  |

|                |                                                                                         |
|----------------|-----------------------------------------------------------------------------------------|
| Product Number | 432027                                                                                  |
| Product Name   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Orange |
| Qty./Pk        | 1 / Pk                                                                                  |
| Qty./Cs        | 1 / Cs                                                                                  |
| Product Number | 432028                                                                                  |
| Product Name   | Corning® CoolBox™ 2XT All-day Cooling and Freezing Workstation, Double Capacity, Pink   |
| Qty./Pk        | 1 / Pk                                                                                  |
| Qty./Cs        | 1 / Cs                                                                                  |

# Corning® Ice Buckets, Round



Color

|  |            |
|--|------------|
|  | Blue       |
|  | Green      |
|  | Orange     |
|  | Pink       |
|  | Purple     |
|  | Red        |
|  | Lime Green |

Capacity

|  |       |
|--|-------|
|  | 2.5 L |
|  | 4 L   |

Products

|                |                                                |
|----------------|------------------------------------------------|
| Product Number | 432122                                         |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Green |
| Qty./Pk        | 1 / Pk                                         |
| Qty./Cs        | 1 / Cs                                         |

|                |                                                     |
|----------------|-----------------------------------------------------|
| Product Number | 432123                                              |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Blue       |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |
| Product Number | 432124                                              |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Red        |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |
| Product Number | 432125                                              |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Orange     |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |
| Product Number | 432126                                              |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Lime Green |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |
| Product Number | 432127                                              |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Pink       |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |
| Product Number | 432128                                              |
| Product Name   | Corning® Ice Bucket with Lid, Round, 4L, Purple     |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |

|                |                                                  |
|----------------|--------------------------------------------------|
| Product Number | 432129                                           |
| Product Name   | Corning® Ice Bucket with Lid, Round, 2.5L, Blue  |
| Qty./Pk        | 1 / Pk                                           |
| Qty./Cs        | 1 / Cs                                           |
| Product Number | 432130                                           |
| Product Name   | Corning® Ice Bucket with Lid, Round, 2.5L, Green |
| Qty./Pk        | 1 / Pk                                           |
| Qty./Cs        | 1 / Cs                                           |
| Product Number | 432131                                           |
| Product Name   | Corning® Ice Bucket with Lid, Round, 2.5L, Red   |
| Qty./Pk        | 1 / Pk                                           |
| Qty./Cs        | 1 / Cs                                           |



# Corning® CoolRack® SV Sample Vial Modules



Cryogenic Storage Format

|                         |
|-------------------------|
| 12 x 10 mL sample vials |
| 12 x 5 mL sample vials  |

Product Type

|               |
|---------------|
| CoolRack SV10 |
| CoolRack SV2  |

## Products

|                |                                                            |
|----------------|------------------------------------------------------------|
| Product Number | 432058                                                     |
| Product Name   | Corning® CoolRack SV2 Holds 12 x 5 mL Injectable Ampules   |
| Qty./Pk        | 1 / Pk                                                     |
| Qty./Cs        | 1 / Cs                                                     |
| Product Number | 432059                                                     |
| Product Name   | Corning® CoolRack SV10 Holds 12 x 10 mL Injectable Ampules |
| Qty./Pk        | 1 / Pk                                                     |
| Qty./Cs        | 1 / Cs                                                     |

# Corning® CoolRack® Tall Modules



**Product Type**

|                    |
|--------------------|
| CoolRack 15 mL     |
| CoolRack 250 mL    |
| CoolRack 250 mL PF |
| CoolRack 50 mL     |
| CoolRack L         |
| CoolRack XT 5 mL   |
| CoolRack XT M24    |

**Cryogenic Storage Format**

|                                                            |
|------------------------------------------------------------|
| 12 x 15 mL centrifuge tubes, light weight insulated module |
| 12 x 5 mL centrifuge tubes, SBS-compatible                 |
| 4 x 50 mL centrifuge tube                                  |
| 9 x 15 mL centrifuge tube                                  |
| One 250 mL centrifuge tube                                 |
| One 250 mL Easy Grip storage bottle                        |
| 24 x 1.5 or 2 mL microcentrifuge tubes, SBS-compatible     |

**Products**

|                |                                                                        |
|----------------|------------------------------------------------------------------------|
| Product Number | 432040                                                                 |
| Product Name   | Corning® CoolRack XT M24, Holds 24 x 1.5 or 2 mL Microcentrifuge Tubes |
| Qty./Pk        | 1 / Pk                                                                 |
| Qty./Cs        | 1 / Cs                                                                 |
| Product Number | 432060                                                                 |
| Product Name   | Corning® CoolRack XT 5mL, Holds 12 x 5mL Microcentrifuge Tubes         |
| Qty./Pk        | 1 / Pk                                                                 |
| Qty./Cs        | 1 / Cs                                                                 |
| Product Number | 432061                                                                 |
| Product Name   | Corning® CoolRack 15mL, Holds 9 x 15mL Conical Centrifuge Tubes        |
| Qty./Pk        | 1 / Pk                                                                 |
| Qty./Cs        | 1 / Cs                                                                 |
| Product Number | 432062                                                                 |
| Product Name   | Corning® CoolRack 50 mL, Holds 4 x 50 mL Conical Centrifuge Tubes      |
| Qty./Pk        | 1 / Pk                                                                 |
| Qty./Cs        | 1 / Cs                                                                 |
| Product Number | 432063                                                                 |
| Product Name   | Corning® CoolRack 250 mL, Holds 1 x 250 mL Conical Centrifuge Tubes    |
| Qty./Pk        | 1 / Pk                                                                 |
| Qty./Cs        | 1 / Cs                                                                 |

|                |                                                                                    |
|----------------|------------------------------------------------------------------------------------|
| Product Number | 432064                                                                             |
| Product Name   | Corning® CoolRack 250 mL, Holds 1 x 250 mL Easy Grip Storage Bottle                |
| Qty./Pk        | 1 / Pk                                                                             |
| Qty./Cs        | 1 / Cs                                                                             |
| Product Number | 432068                                                                             |
| Product Name   | Corning® CoolRack L, Light Weight, Insulated Module for 12 x 15mL Centrifuge Tubes |
| Qty./Pk        | 1 / Pk                                                                             |
| Qty./Cs        | 1 / Cs                                                                             |

# Corning® CoolRack® V Blood Collection Tube Modules



Product Type

|  |               |
|--|---------------|
|  | CoolRack LV   |
|  | CoolRack V13  |
|  | CoolRack V16  |
|  | CoolRack VS13 |

Cryogenic Storage Format

|  |                                                                                                                   |
|--|-------------------------------------------------------------------------------------------------------------------|
|  | 12 x 13 mm or 16 mm diameter tubes, light weight insulated module, thermo-conductive base and insulative exterior |
|  | Nine 13 x 100 mm blood tubes or 5 mL cryogenic vials                                                              |
|  | Nine 13 x 75 mm blood tubes                                                                                       |
|  | Nine 16 x 100 mm blood tubes                                                                                      |

Products

|                |                                                     |
|----------------|-----------------------------------------------------|
| Product Number | 432065                                              |
| Product Name   | Corning® CoolRack VS13, Holds 9x13x75mm Blood Tubes |
| Qty./Pk        | 1 / Pk                                              |
| Qty./Cs        | 1 / Cs                                              |

|                       |                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------|
| <b>Product Number</b> | 432066                                                                                      |
| <b>Product Name</b>   | Corning® CoolRack V16, Holds 9x16x100 mm Blood Tubes                                        |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                      |
| <b>Product Number</b> | 432067                                                                                      |
| <b>Product Name</b>   | Corning® CoolRack V13, Holds 9x13x100 mm Blood Tubes                                        |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                      |
| <b>Product Number</b> | 432069                                                                                      |
| <b>Product Name</b>   | Corning® CoolRack LV, Light Weight, Insulated Module for 12 x 13 mm or 16 mm Diameter Tubes |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                      |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                      |

# Corning® CoolRack® CF Cryogenic Vial and FACS Tube Modules



Cryogenic Storage Format

|                                                                                                           |
|-----------------------------------------------------------------------------------------------------------|
| 24 cryogenic vials or FACS tubes, with "gripping" wells for one-hand vial opening/closing, SBS-compatible |
| 30 cryogenic vials or FACS tubes, with "gripping" wells for one-hand vial opening/closing                 |
| 15 cryogenic vials or FACS tubes                                                                          |
| 45 cryogenic vials or FACS tubes                                                                          |

Product Type

|                   |
|-------------------|
| CoolRack CF15     |
| CoolRack CFT30    |
| CoolRack CF45     |
| CoolRack XT CFT24 |

Products

|                |                                                                |
|----------------|----------------------------------------------------------------|
| Product Number | 432049                                                         |
| Product Name   | Corning® CoolRack CF15, Holds 15 Cryogenic Vials or FACS Tubes |
| Qty./Pk        | 1 / Pk                                                         |
| Qty./Cs        | 1 / Cs                                                         |

|                       |                                                                                                                            |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Number</b> | 432050                                                                                                                     |
| <b>Product Name</b>   | Corning® CoolRack XT CFT24, Holds 24 Cryogenic vial or FACS Tubes, with "Gripping" Wells for One-hand Vial Opening/Closing |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                                     |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                                                     |
| <b>Product Number</b> | 432051                                                                                                                     |
| <b>Product Name</b>   | Corning® CoolRack CF45, Holds 45 Cryogenic Vials or FACS Tubes                                                             |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                                     |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                                                     |
| <b>Product Number</b> | 432052                                                                                                                     |
| <b>Product Name</b>   | Corning® CoolRack CFT30, Holds 30 Cryogenic Vials or FACS Tubes, with "Gripping" Wells for One-hand Vial Opening/Closing   |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                                     |
| <b>Qty./Cs</b>        | 1 / Cs                                                                                                                     |



# Corning® CoolRack® PCR Microplate, Stripwell® or Tube Modules



**Product Type**

|                                   |
|-----------------------------------|
| CoolRack XT PCR96                 |
| CoolRack XT PCR384                |
| CoolRack XT-M-PCR, SBS-Compatible |

**Cryogenic Storage Format**


|                                                                                     |
|-------------------------------------------------------------------------------------|
| 384-well PCR microplate, SBS-compatible                                             |
| 96-well PCR microplate, PCR 12-strip tubes or microcentrifuge tubes, SBS-compatible |
| 12 x 1.5 or 2 mL microcentrifuge tubes or 6-stripwells                              |

**Products**

|                |                                                  |
|----------------|--------------------------------------------------|
| Product Number | 432053                                           |
| Product Name   | Corning® CoolRack XT PCR96, Holds 12 Strip Wells |
| Qty./Pk        | 1 / Pk                                           |
| Qty./Cs        | 1 / Cs                                           |

|                |                                                                                            |
|----------------|--------------------------------------------------------------------------------------------|
| Product Number | 432054                                                                                     |
| Product Name   | Corning® CoolRack XT-M-PCR, Holds 12 x 1.5 or 2 mL Microcentrifuge Tubes and 6 Strip Wells |
| Qty./Pk        | 1 / Pk                                                                                     |
| Qty./Cs        | 1 / Cs                                                                                     |
| Product Number | 432055                                                                                     |
| Product Name   | Corning® CoolRack XT PCR384 for 384-well Plate                                             |
| Qty./Pk        | 1 / Pk                                                                                     |
| Qty./Cs        | 1 / Cs                                                                                     |

# Corning® CoolRack® and CoolSink® Plate and Reservoir Modules



Product Type

|  |                 |
|--|-----------------|
|  | CoolSink LX55   |
|  | CoolSink XT 96F |
|  | CoolSink XT 96U |
|  | CoolRack LX50   |

Cryogenic Storage Format

|  |                                                                                       |
|--|---------------------------------------------------------------------------------------|
|  | One 55 mL reagent reservoir                                                           |
|  | One 6-, 12-, 24-, 48-, 96-, or 384-well flat bottom plate, SBS-compatible             |
|  | One 96-well U-bottom microplate, SBS-compatible                                       |
|  | One Axygen® 50 mL disposable reagent reservoir (RES-V-50, RES-V-50-S, or RES-V-50-SI) |

Products

|                |                                          |
|----------------|------------------------------------------|
| Product Number | 432070                                   |
| Product Name   | CoolSink XT96F, Flat Bottom Plate Module |
| Qty./Pk        | 1 / Pk                                   |
| Qty./Cs        | 1 / Cs                                   |

|                |                                                |
|----------------|------------------------------------------------|
| Product Number | 432071                                         |
| Product Name   | CoolSink XT 96U, 96 Well U-bottom Plate Module |
| Qty./Pk        | 1 / Pk                                         |
| Qty./Cs        | 1 / Cs                                         |
| Product Number | 432072                                         |
| Product Name   | CoolSink LX55, 55mL Reagent Reservoir Module   |
| Qty./Pk        | 1 / Pk                                         |
| Qty./Cs        | 1 / Cs                                         |
| Product Number | 432139                                         |
| Product Name   | CoolRack® LX50 Reservoir Module                |
| Qty./Pk        | 1 / Pk                                         |
| Qty./Cs        | 1 / Cs                                         |

# Corning® CoolBox™ XT Starter Ice-free Coolers



Color

Purple

Cryogenic Storage Format

XT Holder and XT Cooling Core

XT Starter Holder only

Product Type

XT Starter, complete

XT Starter, holder only

## Products

|                |                                                                          |
|----------------|--------------------------------------------------------------------------|
| Product Number | 432014                                                                   |
| Product Name   | Corning® XT Starter Complete, Low profile, Open Platform Ice-free Cooler |
| Qty./Pk        | 1 / Pk                                                                   |
| Qty./Cs        | 1 / Cs                                                                   |
| Product Number | 432015                                                                   |
| Product Name   | Corning® XT Starter Holder Only, Open Platform Ice-free Cooler           |
| Qty./Pk        | 1 / Pk                                                                   |
| Qty./Cs        | 1 / Cs                                                                   |

# Corning® ThermalTray™ Thermo-conductive Platforms



For Use With

|                                              |
|----------------------------------------------|
| For use with 9L ice pan with crushed ice     |
| For use with 9L ice pan with liquid nitrogen |
| Water bath                                   |

Product Type

|                                     |
|-------------------------------------|
| Thermal Tray HP, high-profile       |
| Thermal Tray LP, low profile        |
| Thermal Tray SLP, slim, low profile |

## Products

|                |                                                 |
|----------------|-------------------------------------------------|
| Product Number | 432073                                          |
| Product Name   | ThermalTray SLP Low Profile for Liquid Nitrogen |
| Qty./Pk        | 1 / Pk                                          |
| Qty./Cs        | 1 / Cs                                          |

|                |                                            |
|----------------|--------------------------------------------|
| Product Number | 432074                                     |
| Product Name   | ThermalTray LP, Medium Profile for Ice     |
| Qty./Pk        | 1 / Pk                                     |
| Qty./Cs        | 1 / Cs                                     |
|                |                                            |
| Product Number | 432075                                     |
| Product Name   | ThermalTrayHP, High Profile for Water Bath |
| Qty./Pk        | 1 / Pk                                     |
| Qty./Cs        | 1 / Cs                                     |

# Corning® BioCoat® Poly-D-Lysine/Laminin 100 mm TC-Treated Culture Dishes, 5/Pack, 10/Case

**Product Number 354455**

Corning® BioCoat™ products are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum-free or serum-containing cultures. Corning offers custom coating capabilities.

Shelf life: >= 3 months from date of shipment

Shipping room temperature.

---

## Details

|                   |                       |
|-------------------|-----------------------|
| Qty./Pk           | 5 / Pk                |
| Qty./Cs           | 10 / Cs               |
| Cell Growth Area  | 58.1 cm <sup>2</sup>  |
| Surface Treatment | TC-Treated            |
| Storage           | 2° - 8°C              |
| Surface Coating   | Poly-D-Lysine/Laminin |

---



# Corning® BioCoat® Poly-D-Lysine Culture Dishes



Poly-D-Lysine (PDL) and Poly-L-Lysine (PLL) are synthetic compounds that enhance cell adhesion and protein absorption by altering surface charges on the culture substrate. In addition to promoting cell adhesion, poly-lysine surface treatments support neurite outgrowth and improve the survival of many central nervous system (CNS) primary cells in culture. As PDL and PLL are synthetic molecules, they do not stimulate biological activity in the cells cultured on them, and they do not introduce impurities carried by natural polymers.

**Applications include:**

- Attachment and spreading of a variety of cell types
- Cell differentiation and neurite outgrowth
- Attachment of fastidious transfected cell lines
- Support survival of primary neurons in culture
- Serum-free or reduced-serum culture

**Source**

- PDL, synthetic (MW 75-150 kD)
- PLL, synthetic (MW 30-70 kD)

**Quality Control**

- Tested for ability to promote firm attachment of RCG cells
- Tested and found negative for bacteria and fungi

**Storage**

- 2°C to 8°C under dry conditions.

**Diameter**

|  |        |
|--|--------|
|  | 35 mm  |
|  | 60 mm  |
|  | 100 mm |
|  | 150 mm |

**Qty/Cs**

|  |       |
|--|-------|
|  | 5/Cs  |
|  | 10/Cs |
|  | 20/Cs |

|        |
|--------|
| 40/Cs  |
| 100/Cs |

# Products

|                |                                                                                          |
|----------------|------------------------------------------------------------------------------------------|
| Product Number | 354467                                                                                   |
| Product Name   | Corning® BioCoat® Poly-D-Lysine 35 mm TC-treated Culture Dishes, 20/Pack, 20/Case        |
| Qty./Pk        | 20 / Pk                                                                                  |
| Qty./Cs        | 20 / Cs                                                                                  |
| Product Number | 354468                                                                                   |
| Product Name   | Corning® BioCoat® Poly-D-Lysine 60 mm TC-treated Culture Dishes, 20/Pack, 20/Case        |
| Qty./Pk        | 20 / Pk                                                                                  |
| Qty./Cs        | 20 / Cs                                                                                  |
| Product Number | 354469                                                                                   |
| Product Name   | Corning® BioCoat® Poly-D-Lysine 100 mm TC-treated Culture Dishes, 10Pack, 10/Case        |
| Qty./Pk        | 10 / Pk                                                                                  |
| Qty./Cs        | 10 / Cs                                                                                  |
| Product Number | 354550                                                                                   |
| Product Name   | Corning® BioCoat® Poly-D-Lysine 150 mm TC-treated Gridded Culture Dishes, 5/Pack, 5/Case |
| Qty./Pk        | 5 / Pk                                                                                   |
| Qty./Cs        | 5 / Cs                                                                                   |

|                       |                                                                                    |
|-----------------------|------------------------------------------------------------------------------------|
| <b>Product Number</b> | 356467                                                                             |
| <b>Product Name</b>   | Corning® BioCoat® Poly-D-Lysine 35 mm TC-treated Culture Dishes, 20/Pack, 100/Case |
| <b>Qty./Pk</b>        | 20 / Pk                                                                            |
| <b>Qty./Cs</b>        | 100 / Cs                                                                           |
| <b>Product Number</b> | 356468                                                                             |
| <b>Product Name</b>   | Corning® BioCoat® Poly-D-Lysine 60 mm Culture Dishes, 20/Pack, 100/Case            |
| <b>Qty./Pk</b>        | 20 / Pk                                                                            |
| <b>Qty./Cs</b>        | 100 / Cs                                                                           |
| <b>Product Number</b> | 356469                                                                             |
| <b>Product Name</b>   | Corning® BioCoat® Poly-D-Lysine 100 mm TC-treated Culture Dishes, 10/Pack, 40/Case |
| <b>Qty./Pk</b>        | 10 / Pk                                                                            |
| <b>Qty./Cs</b>        | 40 / Cs                                                                            |

# Corning® BioCoat® Matrigel® Matrix Culture Dishes



Corning Matrigel basement membrane matrix is a solubilized basement membrane preparation extracted from the Engelbreth-Holm-Swarm (EHS) mouse sarcoma, a tumor rich in ECM proteins. Its major component is laminin, followed by collagen IV, heparan sulfate proteoglycans, entactin, and nidogen. Corning Matrigel matrix is effective for the attachment and differentiation of both normal and transformed anchorage-dependent epithelial and other cell types including neurons and oligodendrocytes.

**Applications include:**

- Elicitation of tissue-specific cellular morphology and protein production in epithelial cells
- Differentiation of endothelial, muscle, and neuronal cells
- Development of three-dimensional matrix model systems

**Source**

- EHS mouse tumor

**Formulation**

- Dulbecco's Modified Eagles' Medium with 50 µg/mL gentamycin. Corning Matrigel matrix is compatible with all culture media.

**Quality Control**

- Tested for ability to promote neurite outgrowth from chick dorsal root ganglia in the absence of NGF
- Tested and found negative for bacteria and fungi

**Storage and Stability**

- -20°C. Keep frozen until use.
- Thin layer cultureware stable at 2°C to 8°C.

**Diameter**

|  |        |
|--|--------|
|  | 60 mm  |
|  | 100 mm |

**Qty/Cs**

|  |       |
|--|-------|
|  | 10/Cs |
|  | 20/Cs |

---

**Products**

|                |                                                                                      |
|----------------|--------------------------------------------------------------------------------------|
| Product Number | 354600                                                                               |
| Product Name   | Corning® BioCoat® Matrigel® Matrix Thin-Layer 100 mm Culture Dishes, 5/Pack, 10/Case |
| Qty./Pk        | 5 / Pk                                                                               |
| Qty./Cs        | 10 / Cs                                                                              |
| Product Number | 354601                                                                               |
| Product Name   | Corning® BioCoat® Matrigel® Matrix Thin-Layer 60 mm Culture Dishes, 5/Pack, 20/Case  |
| Qty./Pk        | 5 / Pk                                                                               |
| Qty./Cs        | 20 / Cs                                                                              |

# Falcon® Cell Culture Dishes



- Flat, optically clear polystyrene surfaces for distortion-free microscopic visualization of cells
- Uniform surface chemistry created by vacuum-gas plasma treatment promotes cell attachment
- Improved handling of small dishes with the unique Easy-Grip design
- Nonpyrogenic
- Sterilized by gamma irradiation
- Packaged in peel-open, medical-style bags
- Crystal-grade virgin polystyrene
- Standard Tissue Culture (TC), Corning® Primaria™ surface treatment and non-TC-treated polystyrene available

## Designed for cell culture

- Flat, distortion-free optics
- Lids designed for optimal gas exchange
- Stacking rings allow for easier stacking and handling
- Vacuum-gas plasma treatment permanently and consistently modifies the cell growth surface
- Standard tissue culture surface is hydrophilic and contains a variety of negatively charged functional groups that support cell attachment and spreading
- Corning Primaria tissue culture surface additionally incorporates nitrogen-containing functional groups and has been shown to support improved attachment, spreading, and differentiation of some cell types

## Diameter

|  |        |
|--|--------|
|  | 35 mm  |
|  | 60 mm  |
|  | 100 mm |
|  | 150 mm |

## Dish Type

|  |                           |
|--|---------------------------|
|  | Easy Grip                 |
|  | Standard                  |
|  | Gridded                   |
|  | Center Well Organ Culture |

## Products

|                |                                                                                         |
|----------------|-----------------------------------------------------------------------------------------|
| Product Number | 353001                                                                                  |
| Product Name   | Falcon® 35 mm TC-treated Easy-Grip Style Cell Culture Dish, 20/Pack, 500/Case, Sterile  |
| Qty./Pk        | 20 / Pk                                                                                 |
| Qty./Cs        | 500 / Cs                                                                                |
| Product Number | 353002                                                                                  |
| Product Name   | Falcon® 60 mm TC-treated Cell Culture Dish, 20/Pack, 500/Case, Sterile                  |
| Qty./Pk        | 20 / Pk                                                                                 |
| Qty./Cs        | 500 / Cs                                                                                |
| Product Number | 353003                                                                                  |
| Product Name   | Falcon® 100 mm TC-treated Cell Culture Dish, 20/Pack, 200/Case, Sterile                 |
| Qty./Pk        | 20 / Pk                                                                                 |
| Qty./Cs        | 200 / Cs                                                                                |
| Product Number | 353004                                                                                  |
| Product Name   | Falcon® 60 mm TC-treated Easy-Grip Style Cell Culture Dish, 20/Pack, 500/Case, Sterile  |
| Qty./Pk        | 20 / Pk                                                                                 |
| Qty./Cs        | 500 / Cs                                                                                |
| Product Number | 353025                                                                                  |
| Product Name   | Falcon® 150 mm TC-treated Cell Culture Dish with 20 mm Grid, 10/Pack, 100/Case, Sterile |
| Qty./Pk        | 10 / Pk                                                                                 |
| Qty./Cs        | 100 / Cs                                                                                |

|                |                                                                                     |
|----------------|-------------------------------------------------------------------------------------|
| Product Number | 353037                                                                              |
| Product Name   | Falcon® 60 mm TC-treated Center Well Organ Culture Dish, 20/Pack, 500/Case, Sterile |
| Qty./Pk        | 20 / Pk                                                                             |
| Qty./Cs        | 500 / Cs                                                                            |



# Corning® 245mm Square BioAssay Dishes



Square bioassay dishes are made from polystyrene and are nonpyrogenic. They are packed with lids and are designed with a stacking bead so that they will stack securely without slipping. The dishes are compatible with automated colony picking instruments.

Height

18 mm

Low Profile

Handle

No Handle


Yes

## Products

|                |                                                                           |
|----------------|---------------------------------------------------------------------------|
| Product Number | 431111                                                                    |
| Product Name   | Corning® 245 mm Square BioAssay Dish with Handles, not TC-treated Culture |
| Qty./Pk        | 4 / Pk                                                                    |
| Qty./Cs        | 16 / Cs                                                                   |
| Product Number | 431272                                                                    |
| Product Name   | Corning® 245 mm Square BioAssay Dish without Handles, not TC-treated      |
| Qty./Pk        | 4 / Pk                                                                    |
| Qty./Cs        | 16 / Cs                                                                   |

|                |                                                                                  |
|----------------|----------------------------------------------------------------------------------|
| Product Number | 431301                                                                           |
| Product Name   | Corning® 245 mm Square Low Profile BioAssay Dish without Handles, not TC-treated |
| Qty./Pk        | 5 / Pk                                                                           |
| Qty./Cs        | 20 / Cs                                                                          |

# Corning® Treated Culture Dishes



- Corning® CellBIND® surface is a novel cell culture treatment that increases surface wettability for more even and consistent cell attachment
- Ultra-Low Attachment surface dishes feature a covalently bound hydrogel layer that minimizes cell attachment, protein absorption, and cellular activation
- 10 dishes/bag are available for 100 mm dishes (Cat. No. 430293)
- 245 mm square dishes offer 500 cm<sup>2</sup> growth surface
- Stacking beads aid in handling
- Vents provide consistent gas exchange
- Manufactured from optically clear virgin polystyrene
- Sterile
- Nonpyrogenic

Dish style (mm) = actual growth surface diameters: 35 mm dish = 34.4 mm; 60 mm dish = 52.1 mm; 100 mm dish = 83.8 mm; 150 mm dish = 139.1 mm.

Cat. no. 431110, Square dishes with interior bottom dimensions of 224 x 224 mm.

## Surface Treatment

|                      |
|----------------------|
| Corning CellBIND     |
| TC-treated           |
| Ultra-Low Attachment |

## Dimensions

|                   |
|-------------------|
| 35 x 10 mm        |
| 60 x 15 mm        |
| 60 x 60 mm        |
| 100 x 20 mm       |
| 150 x 25 mm       |
| 245 x 245 x 25 mm |

## Grid

|        |
|--------|
| No     |
| Yes    |
| Qty/Pk |
| 10/Pk  |
| 20/Pk  |
| 4/Pk   |
| 5/Pk   |
| 7/Pk   |

### Products

|                |                                                  |
|----------------|--------------------------------------------------|
| Product Number | 3261                                             |
| Product Name   | Corning® 60 mm Ultra-Low Attachment Culture Dish |
| Qty./Pk        | 5 / Pk                                           |
| Qty./Cs        | 20 / Cs                                          |
| Product Number | 3294                                             |
| Product Name   | Corning® CellBIND® Surface 35 mm Culture Dish    |
| Qty./Pk        | 10 / Pk                                          |
| Qty./Cs        | 210 / Cs                                         |
| Product Number | 3295                                             |
| Product Name   | Corning® CellBIND® Surface 60 mm Culture Dish    |
| Qty./Pk        | 7 / Pk                                           |
| Qty./Cs        | 126 / Cs                                         |
| Product Number | 3296                                             |
| Product Name   | Corning® CellBIND® Surface 100 mm Culture Dish   |
| Qty./Pk        | 5 / Pk                                           |
| Qty./Cs        | 40 / Cs                                          |

|                |                                                       |
|----------------|-------------------------------------------------------|
| Product Number | 430165                                                |
| Product Name   | Corning® 35 mm TC-treated Culture Dish                |
| Qty./Pk        | 20 / Pk                                               |
| Qty./Cs        | 500 / Cs                                              |
| Product Number | 430166                                                |
| Product Name   | Corning® 60 mm TC-treated Culture Dish                |
| Qty./Pk        | 20 / Pk                                               |
| Qty./Cs        | 500 / Cs                                              |
| Product Number | 430167                                                |
| Product Name   | Corning® 100 mm TC-treated Culture Dish               |
| Qty./Pk        | 20 / Pk                                               |
| Qty./Cs        | 500 / Cs                                              |
| Product Number | 430196                                                |
| Product Name   | Corning® 60 mm with 2 mm Grid TC-treated Culture Dish |
| Qty./Pk        | 20 / Pk                                               |
| Qty./Cs        | 500 / Cs                                              |
| Product Number | 430293                                                |
| Product Name   | Corning® 100 mm TC-treated Culture Dishes, 10/Pack    |
| Qty./Pk        | 10 / Pk                                               |
| Qty./Cs        | 480 / Cs                                              |
| Product Number | 430599                                                |
| Product Name   | Corning® 150 mm TC-treated Culture Dish               |
| Qty./Pk        | 5 / Pk                                                |
| Qty./Cs        | 60 / Cs                                               |

# Corning® BioCoat® Collagen I-coated Culture Dishes



Collagen I, found in most tissues and organs, is most plentiful in dermis, tendon, and bone. It is an integral part of the framework that holds cells and tissues together and has been recognized as a useful matrix for improving cell culture. *in vitro* use of collagen can exert effects on the adherence, morphology, growth, migration, and differentiation of a variety of cell types.

**Applications include:**

- Promotion of cell attachment and spreading
- Rapid expansion of cell populations
- Serum-free or reduced serum culture
- Cell adhesion assays
- Improving survival of primary cells in culture

**Source**

- Rat tail tendon

**Quality Control**

- Tested for ability to promote attachment and spreading of HT-1080 human fibrosarcoma cells
- Tested and found negative for bacteria and fungi
- Collagen I purity >90% by SDS-PAGE

**Storage**

- 4°C to 30°C under dry conditions.

**Qty/Cs**

|  |        |
|--|--------|
|  | 5/Cs   |
|  | 10/Cs  |
|  | 20/Cs  |
|  | 40/Cs  |
|  | 100/Cs |

**Diameter**

|  |        |
|--|--------|
|  | 35 mm  |
|  | 60 mm  |
|  | 100 mm |
|  | 150 mm |

Products

|                |                                                                                 |
|----------------|---------------------------------------------------------------------------------|
| Product Number | 354401                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 60 mm TC-treated Culture Dishes, 20/Pack, 20/Case  |
| Qty./Pk        | 20 / Pk                                                                         |
| Qty./Cs        | 20 / Cs                                                                         |
| Product Number | 354450                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 100 mm TC-treated Culture Dishes, 10/Pack, 10/Case |
| Qty./Pk        | 10 / Pk                                                                         |
| Qty./Cs        | 10 / Cs                                                                         |
| Product Number | 354456                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 35 mm TC-treated Culture Dishes, 20/Pack, 20/Case  |
| Qty./Pk        | 20 / Pk                                                                         |
| Qty./Cs        | 20 / Cs                                                                         |
| Product Number | 354551                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 150 mm TC-treated Culture Dishes, 5/Pack, 5/Case   |
| Qty./Pk        | 5 / Pk                                                                          |
| Qty./Cs        | 5 / Cs                                                                          |
| Product Number | 356401                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 60 mm TC-treated Culture Dishes, 20/Pack, 100/Case |
| Qty./Pk        | 20 / Pk                                                                         |
| Qty./Cs        | 100 / Cs                                                                        |

|                |                                                                                 |
|----------------|---------------------------------------------------------------------------------|
| Product Number | 356450                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 100 mm TC-treated Culture Dishes, 10/Pack, 40/Case |
| Qty./Pk        | 10 / Pk                                                                         |
| Qty./Cs        | 40 / Cs                                                                         |
| Product Number | 356456                                                                          |
| Product Name   | Corning® BioCoat® Collagen I 35 mm TC-treated Culture Dishes, 20/Pack, 100/Case |
| Qty./Pk        | 20 / Pk                                                                         |
| Qty./Cs        | 100 / Cs                                                                        |



# Corning® Not TC-Treated Culture Dishes

- Manufactured from optically clear virgin polystyrene
- Not treated for applications where cell attachment is not desired
- Not treated for applications where cell attachment is not desired
- Stacking beads aid in handling
- Vents provide consistent gas exchange
- Sterile
- Sterile
- Nonpyrogenic

Dish style (mm) = actual growth surface diameters: 35 mm dish = 34.4 mm; 60 mm dish = 52.1 mm; 100 mm dish = 83.8 mm; 150 mm dish = 139.1 mm.

Cat. no. 431111, Square dishes with interior bottom dimensions of 224 x 224 mm.

## Dimensions

|             |
|-------------|
| 35 x 10 mm  |
| 60 x 15 mm  |
| 100 x 20 mm |
| 150 x 25 mm |

## Products

|                |                                            |
|----------------|--------------------------------------------|
| Product Number | 430588                                     |
| Product Name   | Corning® 35 mm Not TC-treated Culture Dish |
| Qty./Pk        | 20 / Pk                                    |
| Qty./Cs        | 500 / Cs                                   |

|                       |                                             |
|-----------------------|---------------------------------------------|
| <b>Product Number</b> | 430589                                      |
| <b>Product Name</b>   | Corning® 60 mm Not TC-treated Culture Dish  |
| <b>Qty./Pk</b>        | 20 / Pk                                     |
| <b>Qty./Cs</b>        | 500 / Cs                                    |
| <b>Product Number</b> | 430591                                      |
| <b>Product Name</b>   | Corning® 100 mm Not TC-treated Culture Dish |
| <b>Qty./Pk</b>        | 20 / Pk                                     |
| <b>Qty./Cs</b>        | 500 / Cs                                    |
| <b>Product Number</b> | 430597                                      |
| <b>Product Name</b>   | Corning® 150 mm Not TC-treated Culture Dish |
| <b>Qty./Pk</b>        | 5 / Pk                                      |
| <b>Qty./Cs</b>        | 60 / Cs                                     |

# Corning® Primaria™ Cell Culture Dishes



The Corning Primaria surface features a unique mixture of oxygen-containing (negatively charged) and nitrogen-containing (positively charged) functional groups on the polystyrene surface. The surface supports the growth of cells that can exhibit poor attachment or limited differentiation potential when cultured on traditional TC surfaces, including neuronal, primary, endothelial, and tumor cells. The surface consistency of each lot is confirmed by electron spectroscopy chemical analysis (ESCA).

- Crystal-grade polystyrene modified by proprietary vacuum-gas plasma treatment process
- Stable, permanent surface modification
- Optically clear
- 
- Packaged in red color-coded, peel-open, medical-style packages
- Sterilized by gamma irradiation
- Nonpyrogenic tested to less than 0.1 EU/mL

**Diameter**

|  |        |
|--|--------|
|  | 35 mm  |
|  | 60 mm  |
|  | 100 mm |

## Products

|                |                                                                                         |
|----------------|-----------------------------------------------------------------------------------------|
| Product Number | 353801                                                                                  |
| Product Name   | Corning® Primaria™ 35 mm x 10 mm Easy Grip Style Cell Culture Dish, 20/Sleeve, 200/Case |
| Qty./Pk        | 20 / Pk                                                                                 |
| Qty./Cs        | 200 / Cs                                                                                |

|                       |                                                                                 |
|-----------------------|---------------------------------------------------------------------------------|
| <b>Product Number</b> | 353802                                                                          |
| <b>Product Name</b>   | Corning® Primaria™ 60 mm x 15 mm Standard Cell Culture Dish, 20/Pack, 200/Case  |
| <b>Qty./Pk</b>        | 20 / Pk                                                                         |
| <b>Qty./Cs</b>        | 200 / Cs                                                                        |
| <b>Product Number</b> | 353803                                                                          |
| <b>Product Name</b>   | Corning® Primaria™ 100 mm x 20 mm Standard Cell Culture Dish, 20/Pack, 200/Case |
| <b>Qty./Pk</b>        | 20 / Pk                                                                         |
| <b>Qty./Cs</b>        | 200 / Cs                                                                        |

# Corning® BioCoat® Gelatin Culture Dishes



Corning BioCoat Gelatin cultureware provides an attachment and growth promoting substrate for the culture of a variety of cell types. Gelatin is commonly used in the culture of vascular endothelial cells, muscle, embryonic stem (ES) cells, and F9 teratocarcinoma cells. It is also suitable for promoting adhesion of transfected cell types. Gelatin is a heterogeneous mixture of water-soluble proteins derived through the hydrolysis of Collagen.

**Applications include:**

- Promotion of cell attachment and spreading
- Culture of normal and transfected F9 teratocarcinoma cells for gene expression studies
- Culture of HUVEC for E-Selectin7 expression and VEGF induction8

**Source**

- Gelatin, porcine

**Quality Control**

- Tested for ability to promote proliferation of HUVECs
- Tested and found negative for bacteria and fungi

**Storage**

- 4°C to 30°C under dry conditions.

**Qty/Cs**

|  |       |
|--|-------|
|  | 10/Cs |
|  | 40/Cs |

## Products

|                |                                                                              |
|----------------|------------------------------------------------------------------------------|
| Product Number | 354653                                                                       |
| Product Name   | Corning® BioCoat® Gelatin 100 mm TC-treated Culture Dishes, 10/Pack, 10/Case |
| Qty./Pk        | 10 / Pk                                                                      |
| Qty./Cs        | 10 / Cs                                                                      |

|                |                                                                              |
|----------------|------------------------------------------------------------------------------|
| Product Number | 356653                                                                       |
| Product Name   | Corning® BioCoat® Gelatin 100 mm TC-treated Culture Dishes, 10/Pack, 40/Case |
| Qty./Pk        | 10 / Pk                                                                      |
| Qty./Cs        | 40 / Cs                                                                      |

# Corning® BioCoat® Collagen IV-coated Culture Dishes



Type IV Collagen is a ubiquitous component in basement membranes and provides the major structural support for this matrix. When the Collagen IV meshwork is assembled, it provides a scaffold for the assembly of other basement membrane components through interactions with laminin, entactin/nidogen, and heparan sulfate proteoglycan. Collagen IV is useful as a substrate for growth of epithelial, endothelial, muscle, and nerve cells. Collagen plays a role in the regulation of cell growth, differentiation and adhesion, as well as tissue formation.

**Applications include:**

- Promotion of cell attachment and spreading
- Cell differentiation and neurite outgrowth
- Increased proliferation of PC12 cells
- Studies of effects of collagen IV on cell behavior
- Cell adhesion assays

**Source**

- Engelbreth-Holm-Swarm (EHS) lathrytic mouse tumor

**Quality Control**

- Tested for ability to initiate differentiation (neurite outgrowth) of NG-108 rat glioma/mouse neuroblastoma cells
- Tested and found negative for bacteria and fungi
- Collagen IV purity >90% by SDS-PAGE

**Storage**

- 2°C to 8°C. Do not freeze.

**Diameter**

|  |        |
|--|--------|
|  | 60 mm  |
|  | 100 mm |

---

**Products**

|                |                                                                                 |
|----------------|---------------------------------------------------------------------------------|
| Product Number | 354416                                                                          |
| Product Name   | Corning® BioCoat® Collagen IV 60 mm TC-treated Culture Dishes, 5/Pack, 20/Case  |
| Qty./Pk        | 5 / Pk                                                                          |
| Qty./Cs        | 20 / Cs                                                                         |
| Product Number | 354453                                                                          |
| Product Name   | Corning® BioCoat® Collagen IV 100 mm TC-treated Culture Dishes, 5/Pack, 10/Case |
| Qty./Pk        | 5 / Pk                                                                          |
| Qty./Cs        | 10 / Cs                                                                         |



# Costar® 60 mm Center Well Culture Dish

Product Number 3260



Product is strictly not for human use. This product is not CE marked. For Research Use Only. Not for diagnostic procedures.

- 20mm diameter center well holds 3 mL of medium while the outer well holds 10 mL
- Height is 15mm
- Treated for optimal cell attachment
- Sterilized by gamma radiation and nonpyrogenic
- For Research Use Only. Not for diagnostic procedures.

## Details

|                |          |
|----------------|----------|
| Product Number | 3260     |
| Qty./Pk        | 20 / Pk  |
| Qty./Cs        | 500 / Cs |
| Brand          | Corning® |

# Corning® BioCoat® Laminin Culture Dishes



Laminin (LM), a major component of basement membranes, is a multifunctional glycoprotein that is used as a substrate to culture and maintain differentiated function of a wide variety of cells. Laminin has been shown in culture to stimulate neurite outgrowth, promote cell attachment, chemotaxis, cell differentiation, and neuronal survival.

**Applications include:**

- Promotion of cell attachment and spreading
- Induction of cell differentiation and neurite outgrowth
- Increases proliferation of myoblasts
- Studies of effects of laminin on cell behavior
- Cell adhesion assays

**Source**

- Engelbreth-Holm-Swarm (EHS) mouse tumor

**Quality Control**

- Tested for ability to initiate differentiation (neurite outgrowth) of NG-108 rat glioma/mouse neuroblastoma cells
- Tested and found negative for bacteria and fungi
- Laminin purity >90% by SDS-PAGE (contains entactin)

**Storage**

- 2°C to 8°C. Do not freeze.

**Diameter**

|        |
|--------|
| 60 mm  |
| 100 mm |

## Products

|                |                                                                            |
|----------------|----------------------------------------------------------------------------|
| Product Number | 354405                                                                     |
| Product Name   | Corning® BioCoat® Laminin 60 mm TC-treated Culture Dishes, 5/Pack, 20/Case |
| Qty./Pk        | 5 / Pk                                                                     |
| Qty./Cs        | 20 / Cs                                                                    |

|                |                                                                             |
|----------------|-----------------------------------------------------------------------------|
| Product Number | 354452                                                                      |
| Product Name   | Corning® BioCoat® Laminin 100 mm TC-treated Culture Dishes, 5/Pack, 10/Case |
| Qty./Pk        | 5 / Pk                                                                      |
| Qty./Cs        | 10 / Cs                                                                     |

# Disposable Anti-static Spatulas



- Corning spatulas are designed to save time and to provide contamination-free samples.
- Eliminates the recycling and resterilizing necessary with reusable spatulas.
- Available in five different configurations.
- Microspatulas available in two configurations.
- Individually packaged.
- Antistatic.
- Sterile.
- RNase-/DNase-free.
- Nonpyrogenic.

## Spatula Type

|                                 |
|---------------------------------|
| Spatula, Flat End/Spoon         |
| Microspatula, Rounded End/Scoop |
| Microspatula, Tapered End/Scoop |
| Spatula, Round End/Spoon        |
| Spatula, Small Spoon/Spoon      |
| Spatula, Tapered Blade/Spoon    |
| Spatula, V-Scoop/Spoon          |

## Products

|                |                                                    |
|----------------|----------------------------------------------------|
| Product Number | 3003                                               |
| Product Name   | Corning® Spatula with Tapered Blade/Spoon, Sterile |
| Qty./Pk        | 1 / Pk                                             |
| Qty./Cs        | 100 / Cs                                           |

|                |                                                                  |
|----------------|------------------------------------------------------------------|
| Product Number | 3004                                                             |
| Product Name   | Corning® Spatula with Narrow Spoon/Spoon, Sterile                |
| Qty./Pk        | 1 / Pk                                                           |
| Qty./Cs        | 100 / Cs                                                         |
| Product Number | 3005                                                             |
| Product Name   | Corning® Spatula with Round End Blade/Spoon, Sterile             |
| Qty./Pk        | 1 / Pk                                                           |
| Qty./Cs        | 100 / Cs                                                         |
| Product Number | 3006                                                             |
| Product Name   | Corning® Spatula with V-Scoop/Spoon, Sterile                     |
| Qty./Pk        | 1 / Pk                                                           |
| Qty./Cs        | 100 / Cs                                                         |
| Product Number | 3007                                                             |
| Product Name   | Corning® Spatula with Flat End Blade/Spoon, Sterile              |
| Qty./Pk        | 1 / Pk                                                           |
| Qty./Cs        | 100 / Cs                                                         |
| Product Number | 3012                                                             |
| Product Name   | Corning® Microspatula with V-Shaped Scoop/Tapered Blade, Sterile |
| Qty./Pk        | 1 / Pk                                                           |
| Qty./Cs        | 50 / Cs                                                          |

|                |                                                                  |
|----------------|------------------------------------------------------------------|
| Product Number | 3013                                                             |
| Product Name   | Corning® Microspatula with V-Shaped Scoop/rounded Blade, Sterile |
| Qty./Pk        | 1 / Pk                                                           |
| Qty./Cs        | 50 / Cs                                                          |

# Corning® Gosselin™ Spatulas



- Polystyrene (PS)
- Packed in tamper-evident resealable bag

Length

|  |       |
|--|-------|
|  | 63 mm |
|  | 85 mm |

## Products

|                |                                                                               |
|----------------|-------------------------------------------------------------------------------|
| Product Number | SPA6-01                                                                       |
| Product Name   | Corning® Gosselin™ Spatula, Length 63 mm, White PS, Sterile, 1/Bag, 2000/Case |
| Qty./Pk        | 2000 / Pk                                                                     |
| Qty./Cs        | 2000 / Cs                                                                     |
| Product Number | SPA8-01                                                                       |
| Product Name   | Corning® Gosselin™ Spatula, Length 85 mm, White PS, Sterile, 1/Bag, 2000/Case |
| Qty./Pk        | 2000 / Pk                                                                     |
| Qty./Cs        | 2000 / Cs                                                                     |

# Corning® Temperature Probe Support Clamp

Product Number 440141



Corning offers a 7/16 (16.5cm) support clamp for Corning Cat. No. 6795PR and 400188 temperature probes. Support has two openings, 1/3 (9mm) and 0.15 (3.75mm), for holding probes. It is designed to be used with a boss head clamp (Corning Cat. No. 440140) or standard right angle clamp.

## Details

|                |          |
|----------------|----------|
| Product Number | 440141   |
| Qty./Pk        | 1 / Pk   |
| Qty./Cs        | 1 / Cs   |
| Brand          | Corning® |



# Corning® 18 x 5/16“ Stainless Steel Support Rod

Product Number 440129



Corning offers a stainless steel, two piece, 18 (0.8cm) diameter that consists of two 9 rods that screw together. It screws into the base of Corning hot plate and stirrer models PC-400/410/420, PC-400D/410D/420D and PC-600D/610D/620D.

## Details

|                |          |
|----------------|----------|
| Product Number | 440129   |
| Qty./Pk        | 1 / Pk   |
| Qty./Cs        | 1 / Cs   |
| Brand          | Corning® |

# Corning® 5/16 x 12 Inch Stainless Steel Support Rod

**Product Number 400084**

These 5/16 x12 inch stainless steel support rods screw into the base of Corning® PC-200/210/220 series hot plate/stirrers.

---

## Details

|                |          |
|----------------|----------|
| Product Number | 400084   |
| Qty./Pk        | 1 / Pk   |
| Qty./Cs        | 1 / Cs   |
| Brand          | Corning® |

---

# Corning® Digital Display Hot Plate and Stirring Hot Plate Accessory Kit

**Product Number 6795KIT**

This accessory kit is designed for all Corning digital display hot plates and stirring hot plates (Models PC-400D, PC-420D, PC-600D, PC-620D).

**This Corning® kit Cat. No. 6795KIT includes:**

- Cat. No. 6795PR, external temperature controller
- Cat. No. 440129, vertical support rod. Supplied as two 9 inch (22.9cm) support rods that can be screwed together (18 inches)
- Cat. No. 6970SR, stir bar retriever

**Available for purchase separately and not included in the kit:**

- Cat. No. 440140, boss head clamp
- Cat. No. 440141, holding rod

---

## Details

|                |          |
|----------------|----------|
| Product Number | 6795KIT  |
| Qty./Pk        | 1 / Pk   |
| Qty./Cs        | 1 / Cs   |
| Brand          | Corning® |

---

# Corning® Boss Head Support Clamp

Product Number 440140



Corning offers a 2 diameter support rods. It is designed to be used with the Corning Temperature probe support clamp (Cat. No. 440141).

Note: Bosshead clamps to support rod; holding rod connects to bosshead. Temperature controller slides into holding rod.

## Details

|                |          |
|----------------|----------|
| Product Number | 440140   |
| Qty./Pk        | 1 / Pk   |
| Qty./Cs        | 1 / Cs   |
| Brand          | Corning® |

# PYREX® VISTA™ Single Metric Volumetric and Graduated Cylinders, Plastic Hexagonal Base, Blue Graduations, To Contain



These cylinders are calibrated “to contain”. They are supplied in two parts a strong, accurate glass graduate and a detachable, sturdy plastic base which absorbs impacts and reduces breakage. Tolerances are  $\pm 5\%$  of total at any point. By removing the base, the graduate can be conveniently stacked. Bumper guards are supplied with 25 mL through 100 mL inclusive. The 10 mL has a funnel top for careful pouring.

Capacity

|       |
|-------|
| 10 mL |
|-------|

Products

|                |                                                                                                                    |
|----------------|--------------------------------------------------------------------------------------------------------------------|
| Product Number | 70075-10                                                                                                           |
| Product Name   | PYREX® VISTA™ 10 mL Single Metric Volumetric and Graduated Cylinders, TC, Plastic Hexagonal Base, Blue Graduations |
| Qty./Pk        | 12 / Pk                                                                                                            |
| Qty./Cs        | 12 / Cs                                                                                                            |

# Corning® Single Metric Scale, Graduated Cylinder, Polypropylene, with Funnel Top, To Contain



Translucent polypropylene graduated cylinders are manufactured with raised graduations at the primary scale points. Cylinders have a hexagonal base and are calibrated “To Contain”. Tolerances according to ISO 6706. Exposure to temperatures above 60°C may affect accuracy.

**Capacity**

|  |        |
|--|--------|
|  | 10 mL  |
|  | 25 mL  |
|  | 50 mL  |
|  | 100 mL |
|  | 250 mL |
|  | 500 mL |
|  | 1 L    |
|  | 2 L    |

Products

|                |                                                                                                             |
|----------------|-------------------------------------------------------------------------------------------------------------|
| Product Number | 3022P-10                                                                                                    |
| Product Name   | Corning® Single Metric Scale, 10 mL Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top  |
| Qty./Pk        | 1 / Pk                                                                                                      |
| Qty./Cs        | 12 / Cs                                                                                                     |
| Product Number | 3022P-100                                                                                                   |
| Product Name   | Corning® Single Metric Scale, 100 mL Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top |
| Qty./Pk        | 1 / Pk                                                                                                      |
| Qty./Cs        | 12 / Cs                                                                                                     |
| Product Number | 3022P-1L                                                                                                    |
| Product Name   | Corning® Single Metric Scale, 1L Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top     |
| Qty./Pk        | 1 / Pk                                                                                                      |
| Qty./Cs        | 6 / Cs                                                                                                      |
| Product Number | 3022P-25                                                                                                    |
| Product Name   | Corning® Single Metric Scale, 25 mL Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top  |
| Qty./Pk        | 1 / Pk                                                                                                      |
| Qty./Cs        | 12 / Cs                                                                                                     |
| Product Number | 3022P-250                                                                                                   |
| Product Name   | Corning® Single Metric Scale, 250 mL Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top |
| Qty./Pk        | 1 / Pk                                                                                                      |
| Qty./Cs        | 6 / Cs                                                                                                      |



|                       |                                                                                                             |
|-----------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Product Number</b> | 3022P-2L                                                                                                    |
| <b>Product Name</b>   | Corning® Single Metric Scale, 2L Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top     |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                      |
| <b>Qty./Cs</b>        | 3 / Cs                                                                                                      |
| <b>Product Number</b> | 3022P-50                                                                                                    |
| <b>Product Name</b>   | Corning® Single Metric Scale, 50 mL Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top  |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                      |
| <b>Qty./Cs</b>        | 12 / Cs                                                                                                     |
| <b>Product Number</b> | 3022P-500                                                                                                   |
| <b>Product Name</b>   | Corning® Single Metric Scale, 500 mL Reusable Plastic Graduated Cylinder, Polypropylene, TC with Funnel Top |
| <b>Qty./Pk</b>        | 1 / Pk                                                                                                      |
| <b>Qty./Cs</b>        | 6 / Cs                                                                                                      |

Corning® 100 mL Polystyrene Graduated Cylinder, Sterile

Product Number 430182



These 100 mL Corning® optically clear polystyrene graduated cylinders are designed for aseptic and accurate dispensing of sterile culture media or other biological fluids. A polyethylene dust cover is included.

Details

|             |                                                                                                               |
|-------------|---------------------------------------------------------------------------------------------------------------|
| Qty./Pk     | 1 / Pk                                                                                                        |
| Qty./Cs     | 50 / Cs                                                                                                       |
| Capacity    | 100 mL                                                                                                        |
| Application | Environmental and Water Testing Food and Beverage<br>Testing Pharmaceutical Quality Control and Life Sciences |

**По вопросам продаж и поддержки обращайтесь:**

Алматы (7273)495-231  
Ангарск (3955)60-70-56  
Архангельск (8182)63-90-72  
Астрахань (8512)99-46-04  
Барнаул (3852)73-04-60  
Белгород (4722)40-23-64  
Благовещенск (4162)22-76-07  
Брянск (4832)59-03-52  
Владивосток (423)249-28-31  
Владикавказ (8672)28-90-48  
Владимир (4922)49-43-18  
Волгоград (844)278-03-48  
Вологда (8172)26-41-59  
Воронеж (473)204-51-73  
Екатеринбург (343)384-55-89  
Иваново (4932)77-34-06  
Ижевск (3412)26-03-58  
Иркутск (395)279-98-46  
Казань (843)206-01-48  
Россия +7(495)268-04-70

Калининград (4012)72-03-81  
Калуга (4842)92-23-67  
Кемерово (3842)65-04-62  
Киров (8332)68-02-04  
Коломна (4966)23-41-49  
Кострома (4942)77-07-48  
Краснодар (861)203-40-90  
Красноярск (391)204-63-61  
Курск (4712)77-13-04  
Курган (3522)50-90-47  
Липецк (4742)52-20-81  
Магнитогорск (3519)55-03-13  
Москва (495)268-04-70  
Мурманск (8152)59-64-93  
Набережные Челны (8552)20-53-41  
Нижний Новгород (831)429-08-12  
Новокузнецк (3843)20-46-81  
Ноябрьск (3496)41-32-12  
Новосибирск (383)227-86-73  
Киргизия +996(312)-96-26-47

Омск (3812)21-46-40  
Орел (4862)44-53-42  
Оренбург (3532)37-68-04  
Пенза (8412)22-31-16  
Петрозаводск (8142)55-98-37  
Псков (8112)59-10-37  
Пермь (342)205-81-47  
Ростов-на-Дону (863)308-18-15  
Рязань (4912)46-61-64  
Самара (846)206-03-16  
Саранск (8342)22-96-24  
Санкт-Петербург (812)309-46-40  
Саратов (845)249-38-78  
Севастополь (8692)22-31-93  
Симферополь (3652)67-13-56  
Смоленск (4812)29-41-54  
Сочи (862)225-72-31  
Ставрополь (8652)20-65-13  
Сургут (3462)77-98-35  
Казахстан +7(7172)727-132

Сыктывкар (8212)25-95-17  
Тамбов (4752)50-40-97  
Тверь (4822)63-31-35  
Тольятти (8482)63-91-07  
Томск (3822)98-41-53  
Тула (4872)33-79-87  
Тюмень (3452)66-21-18  
Ульяновск (8422)24-23-59  
Улан-Удэ (3012)59-97-51  
Уфа (347)229-48-12  
Хабаровск (4212)92-98-04  
Чебоксары (8352)28-53-07  
Челябинск (351)202-03-61  
Череповец (8202)49-02-64  
Чита (3022)38-34-83  
Якутск (4112)23-90-97  
Ярославль (4852)69-52-93